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**ENGLISH****INTENDED USE****For in vitro diagnostic use**

Owrens buffer GHI150 is intended for dilution of samples in connection with laboratory diagnostic analysis valid for the coagulation system, mainly when analysing Prothrombin Complex Activity (Owren´s PT-activity). Samples can be citrated plasma, capillary blood and citrated blood.

Owren´s buffer GHI150 may also be used when determining individual coagulation factors such as VIII, IX, VII and X and also when analysing Fibrinogen, Antithrombin and D-dimer.

**PRODUCT DESCRIPTION**

GHI150 is a salt solution with 11mM tri-Sodium Citrate that is buffered to pH 7,3 with 4mM Diethyle Barbital Buffer and with additives of 132 mM NaCl as well as 4 mM Sodium Azid as a conservative.

GHI150 contains 6 x 100 mL Buffer, ready to use.

**PRECAUTIONS**

The product contains Sodium Azid to prevent bacterial growth. Do not empty into drains. Avoid contact with skin and eyes. For more information refer to Material Safety Data Sheet

**REKONSTITUTION**

The product is ready for use.

**STORAGE CONDITIONS AND STABILITY**

Unopened Owren´s Buffer is stable until the expiration date shown on the vial when stored at 15-25°C.

Opened Owren´s Buffer is stable until the expiration date shown on the vial when stored at 15-25°C provided no contamination occurs.

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