Calibrator Plasma Art.No: MRX1201**ENGLISH****INTENDED USE**

For the calibration of MRX Antithrombin liquid reagent.

SUMMARY AND PRINCIPLE

The calibrator plasma is prepared from citrated human plasma with specified level of antithrombin.
The calibrator plasma may be used to generate a standard reference curve for antithrombin determination.

PRODUCT DESCRIPTION

The calibration plasma MRX1201 consists of 1 x 1 mL lyophilized citrated human plasma with specified level of antithrombin.
The calibration value are specified in the lot specific certificate and on the vial label.
The reported value is determined over multiple runs using specific lot of MRX1200 and against SSC/ISTH Secondary Coagulation Standard. SSC/ISTH secondary coagulation standard is traceable to WHO international standards.

PRECAUTIONS*Only For in vitro Diagnostic Use*

The product contains material of human origin. The plasma used in the production is tested and found nonreactive for Hepatitis B surface antigen (HBsAg), Anti-HCV and HIV antibodies. No test can however completely exclude the presence of infected material and the product should be treated as potentially infectious. Waste is disposed of according to local regulations. Wear appropriate clothing. Avoid contact with skin and eyes.

PREPARATION

Allow the vial of calibrator to equilibrate at 15-25°C for 10-15 minutes before opening and reconstitution. *Dissolve the content of each vial with 1.0 mL of CLSI CLRW type water or equivalent¹.* The water temperature should be 15 to 25°C.
Replace the stopper and swirl gently.

Keep the reconstituted calibrator at 15-25°C for 15-30 minutes and mix properly before use. Make sure of the complete reconstitution of the product before use.

STORAGE CONDITIONS AND STABILITY

Unopened calibrators are stable until the expiration date shown on the vial when stored at 2-8°C.

Stability after reconstitution: 8 hours at 15-25°C in closed original vial.

INSTRUMENT/TEST PROCEDURE

Standard curve:

Three points is recommended. (Approximately 0,25 – 0,50 – 1,00 IU/mL)
It is recommended that the analytical equipment automatically dilute the standard curve. Refer to the instruction manual of the manufacturer.
The user must complete a standard curve for each new lot of reagents and / or if control plasma falls outside the assigned limits.
For detailed instructions on performing the tests see instrument manual.

QUALITY CONTROL

In accordance with good laboratory practice it is necessary to run controls to ensure accuracy and reproducibility of the results. It is recommended to use two or three different levels of MediRox controls.
3-Level controls MRX170-MRX183 or
2-Level controls GHI162-GHI170

LIMITATIONS

This product is designed for the calibration of coagulation assays. The calibration plasma is subjected to the limitations of the assay system. Deviations may indicate possible problems with one or more components in the test system.

EXTERNAL PRODUCTS THAT ARE NOT INCLUDED IN THE CALIBRATION KIT MRX 1201

MRX1200 MRX Antithrombin liquid reagent
MRX170-MRX183 3-Level Controls
GHI162-GHI170 2-Level Controls

Revision: 2016-03-14

Changes: New layout, logotype and headings.

Clarification of reconstitution and storage and precautions.

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