PT buffer, polybrene Art.No: MRX152

SVENSKA

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For In Vitro Diagnostic Use

INTENDED USE

PT buffer MRX152 is intended for dilution of samples when analysing Prothrombin complex Activity (Owren's PT-activity). Samples can be citrated plasma, capillary blood and citrated blood. PT Buffer MRX52 contains 5mg/L Polybrene and therefore neutralizes Heparin.

The Buffer can also be used for dilution of samples when determining individual coagulation factors, such as VII, VIII, IX and X and for the determining of Fibrinogen, Antithrombin and D-Dimer. However, this can be reagent dependent and each user should verify the proper functioning of each case.

PRODUCT DESCRIPTION

MRX52 is a Sodium Chloride solution containing 11mM tri-Sodium Citrate, 5mg/L Polybrene, 19mM Hepes, 132 mM NaCl, and 4 mM Sodium Azide as a conservative. pH 7,3.

MRX52 contains 10 x 10 mL Buffer, ready to use.

PRECAUTIONS

The product contains Sodium Azide (<0,1%) to prevent bacterial growth. Do not empty into drains. Avoid contact with skin and eyes.

For more information refer to Material Safety Data Sheet

REKONSTITUTION

The product is ready for use

STORAGE CONDITIONS AND STABILITY

Unopened PT Buffer is stable until the expiration date shown on the vial when stored at 15-25°C.

Opened PT Buffer is stable until the expiration date shown on the vial when stored at 15-25°C provided no contamination occurs.

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