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**ENGLISH****Only For In Vitro Diagnostic Use****INTENDED USE**

PT Buffer MRX150 is intended for dilution of samples when analysing Prothrombin Complex Activity (Owren's PT activity). Samples can be citrated plasma, capillary blood and citrated blood.

The Buffer can also be used for dilution of samples when determining individual coagulation factors, such as VII, VIII, IX and X and for the determining of Fibrinogen, Antithrombin and D-Dimer. However this can be reagent dependent and each user should verify the proper functioning of each case.

**PRODUCT DESCRIPTION**

MRX150 is a Sodium Chloride solution with 11mM tri-Sodium Citrate, 19mM Hepes, 132 mM NaCl and 4 mM Sodium Azid as a conservative. pH 7,3.  
MRX150 contains 10 X 10 mL Buffer, ready to use.

**PRECAUTIONS**

The product contains Sodium Azid (<0,1%) to prevent bacterial growth. Do not empty into drains.

Avoid contact with skin and eyes.

For more information refer to Material Safety Data Sheet

**RECONSTITUTION**

The product is ready for use.

**STORAGE CONDITIONS AND STABILITY**

Unopened PT Buffer is stable until the expiration date shown on the vial when stored at 15-25°C.

Opened PT Buffer is stable until the expiration date shown on the vial when stored at 15-25°C provided no contamination occurs.

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