PT Control INR >5 Art.No: MRX168**SVENSKA**

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For In Vitro Diagnostic use.

INTENDED USE

MRX168 Coagulation Control is intended for use as quality control in the analysis of prothrombin (PT) in the higher measurement range (4-6 INR). The product can be used for the measurement of both PT Owren and PT Quick reagent.

PRODUCT DESCRIPTION

Package: 10 x 1 mL

MRX168 is a freeze-dried preparation of processed citrate anti coagulated human plasma with PT Owren level in the abnormal range.

PRECAUTIONS

The calibrator plasma is of human origin and should be handled as potentially infectious. Each unit human plasma used to manufacture this control has been tested and found non-reactive for Hepatitis B surface Antigen (HBsAg), Anti-HCV and HIV-1/HIV-2 antibodies.

All human source material should be considered as potentially infectious and should be handled with the same precautions as for patient specimens.

Dispose of all waste materials according to the local regulations.

Wear suitable clothing for protection. Avoid contact with skin and eyes.

RECONSTITUTION

Reconstitute each vial with 1.0 mL of distilled or deionized water. The water temperature should be 15-25°C. Replace the stopper and allow the control to stand for 15-30 minutes, swirling occasionally.

Before use invert the vial several times and control visually that all material is completely dissolved.

PROCEDURE

MRX168 should after reconstitution be handled in the same manner as patient specimens and run in accordance with the instructions accompanying the instrument and reagent being used.

STABILITY AND STORAGE

Unopened reagent is stable until the expiration date shown on the label when stored at 2-8°C.

Stability after reconstitution: 5 hours at 15-25°C in the original vial.

EXPECTED VALUES

For lotspecific results see certificate.

PT reagents from different manufacturers may have wide variation at high INR. Both the reference range and expected variation must therefore be assigned and verified by each user. Referenced values below are only to see as a guideline for the expected result.

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