Eximius Control Plus Haemostasis Controls Art.No MRX181

ENGLISH

INTENDED USE

Eximius Control PLUS Level1 is intended for use as quality control plasma for coagulation assays in the normal range.

SUMMARY AND PRINCIPLE

The control is intended for the assessment of precision and accuracy of coagulation analyses such as PT (Quick and Owren method), APTT, Fibrinogen, Antithrombin, Thrombin Time, D-dimer, Plasminogen, Protein C, Protein S and Factors. All analytes are within the normal range except for D-dimer which is slightly increased.

PRODUCT DESCRIPTION

MRX181Eximius Control PLUS kit consists of:

- MRX181 Normal Control Level 1 10 x 1mL lyophilized plasma

This product is prepared from human blood plasma, with added purified biochemicals and preservatives. The controls are provided in lyophilized form for increased stability.

PRECAUTIONS

Biological source material. Handle as potentially infectious. Each human donor unit used to manufacture this control was tested by CE accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), Anti-HCV and HIV-antibodies.

This product is for In Vitro Diagnostic Use.

RECONSTITUTION

Reconstitute each vial with 1.0 mL of distilled or deionized water. The water temperature should be 15-25°C. Replace the stopper and allow the control to stand for 15-30 minutes, swirling occasionally.

Before use gently invert the vial several times and verify visually that all material is completely dissolved.

STORAGE AND STABILITY

Unopened control is stable until the expiration date shown on the vial when stored at 2-8°C. Stability after reconstitution at 2-25°C in sealed original vial:

- Fibrinogen, APTT, Protein S, Protein C, Thrombin Time, Antithrombin, Plasminogen, and D-dimer 54 hours
- PT (Owren and Quick) 24 hours

Factors should be analysed within 1 hour after reconstitution.

To freeze a reconstituted control is common practice, but due to different freezing and thawing processes and other environmental factors that may affect the product, it is up to each user to verify functionality and stability of the product after such procedure.

PROCEDURE

After reconstitution the control should be handled in the same manner as a fresh citrated plasma and run in accordance with the instructions accompanying the instrument, kit, or reagent being used.

LIMITATIONS

This control is designed as control plasma for the monitoring of coagulation assays. The product is not intended for use as a standard/calibrator. Deviations may indicate possible problems with one or more components in the test system.

ASSIGNED VALUES

Mean values and ± 2 SD ranges printed in this insert are determined by User Groups with specific reagents and instruments, in combination with mathematical correction algorithms. The reported values should only be considered as indicative and each laboratory should establish own means and reference ranges.

Note that results may differ considerably depending on laboratory techniques, instrument settings and reagent depending factors.

The control values should be re-established whenever any of the system variables is modified.

PERFORMANCE CHARACTERISTICS

Precision: Within run precision was assessed over multiple runs using specific lots of reagents and control. The coefficient of variation obtained in this study was:

$\leq 1\%$ for Thrombin Time, PT Owren, PT Quick, Free Protein S, Plasminogen; $\leq 2\%$ for, APTT, Protein C, Antithrombin; $\leq 3\%$ for D-dimer $\leq 5\%$ for Fibrinogen.

This product is a freeze-dried product manufactured under rigid quality control standards. To obtain consistent vial-to-vial assay values, the control requires proper storage and handling as described.

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