Calcium chloride Art.No: GHI155, GHI155-2, GHI155-4, GHI155-5, GHI155-10**SVENSKA****Kalciumklorid**

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For In Vitro Diagnostic Use

INTENDED USE

Reagent intended to be used in laboratory diagnostic analyses valid for the coagulationsystem. The product is designed to increase calcium ion activity, when analysing tissue-factor induced coagulation in blood-plasma from citrate-anticoagulated blood, when determining the activity of prothrombin-time (PT).

GHI 155 can also be used when determining Activated Partial Thromboplastin Time (APTT) and specific coagulation factors, such as factor VIII, IX, VII and X.

PRODUCT DESCRIPTION

GHI 155 contains: 25 mM CaCl₂ ready for use solution. The product contains 4 mM NaN₃ as preservative.

GHI155 contains: 6 x 100 mL ready for solution.

GHI155-2 contains 10 x 2 mL ready for use solution.

GHI155-4 contains: 10 x 4 mL ready for use solution.

GHI155-5 contains: 10 x 5 mL ready for use solution.

GHI155-10 contains: 10 x 10 mL ready for use solution.

PRECAUTIONS

The product contains Sodium Azid (<0,1%) to prevent bacterial growth. Do not empty into drains. Avoid contact with skin and eyes.

For more information refer to Material Safety Data Sheet

RECONSTITUTION

The product is ready for use.

For method specification, reference is made to the instrument specific application for each analysis.

STABILITY AND STORAGE

Unopened vials are stable until the expiration date printed on the label when stored at 15-25°C. Opened vials are stable until the expiration date printed on the label if no contamination occurs. Do not store in direct sunlight.

Revision: 2016-04-05

Changes: New logotype and text format All GHI155 products combined into one insert. GHI155-4 and GHI155-10 new. Text about precautions updated.

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