

## MediRox, a coagulation partner to trust



### Company introduction

MediRox has produced instruments, diagnostic reagents and quality control material used in the areas of visco-elasticity and coagulation diagnostics since 1998. Our customers are located within laboratory diagnostics in healthcare and industrial applications.

The company is formed out of a long tradition of scientific coagulation and viscosity research, and is a merger between the Global Hemostasis Institute -a research company originating from Linköping's University Hospital- and the manufactory plant of the well known Chromogenix AB.

MediRox is ISO9000 certified, and all products are CE marked. The company is a public company noted in Stockholm.

### What MediRox can offer:

#### *-National Distribution.*

MediRox is constantly searching for dedicated distribution partners for our product portfolio.

#### *-OEM/ Private labelling/Bulk*

MediRox has a long tradition of OEM partnerships and bulk supply to well known companies including Siemens/DADE, Trinity, IL (ILS Scandinavia), Teco, Helena Lab and others.

*For Products listing in short refer to reverse side or visit [www.medirox.se](http://www.medirox.se)*

## Products listing

### *Coagulation Multi Controls*

MediRox is the first to manage to include all common parameters in one lyophilized coagulation control, including D-Dimer and protein S and C, a product which can save the lab from buying a number of expensive controls.

- **Control plasma.** (Normal and abnormal level)  
Certified levels of PT, APTT, Fibrinogen and Antithrombin.
- **Multi Control plasma.** (Normal, abnormal and extreme abnormal level)  
Certified levels of PT, APTT, Fibrinogen, Antithrombin and D-Dimer.
- **Multi Control plasma PLUS.** (Normal and abnormal level)  
Certified levels of PT, APTT, Fibrinogen, Antithrombin, D-Dimer and proteins S and C

### *Coagulation Reagents*

- **PT Prothrombin Time.**  
PT Owren including buffer, diluents and calcium.  
PT Quick.(ISI=1 and ISI=2) including diluents.
- **The world's first multi platform D-Dimer**  
D-Dimer is a common coagulation analysis on both clinical chemistry and coagulation instruments. MediRox is the first to offer a Multi platform D-Dimers suitable to most instruments on the market. (350-800nm light and 550-650 scattering systems).
- **AT-III LR.** Liquid reagent ready to use.  
Liquid and ready to use Antithrombin reagent.
- **APTT LR** Liquid reagent ready to use.  
Lupus insensitive APTT test for monitoring heparin treated patients and screening of intrinsic factor deficiencies.

### *Instruments*

#### ***-New instrument for advanced coagulation analysis***

Thromboelastography is a common analysis typically used in advanced coagulation investigations, thorax surgery and IR units. MediRox new **ReoRox G<sub>2</sub>** is a modern update to the technology used in ROTEM® or TEG® analysers taking Thromboelastography to new levels.

#### Available systems

- **ReoRox G<sub>2</sub>-4.** A 4-channel device primarily intended for advanced coagulation investigations and as a bed side tool during surgery.
- **ReoRox G<sub>2</sub>-1.** A small one-channel stand alone device intended for general visco-elasticity measurements. Suitable both for health care applications such as whole blood, plasma and other body fluids- as well as general viscosity measurements within the industry.